# FITNESS STAR A PERSONAL FITNESS ASSISTANT EXPERT SYSTEM

# **INDEX**

Topic	Page number
Abstract	3
Features of the Fitness Star Expert System	3
Rules and Descriptions	4
Usage Manual	6
Sample runs (Run #1, #2 and #3)	7
Test cases	10

**Note:** Suitable links have been provided for additional information wherever necessary in RULES and DESCRIPTION.

#### Abstract:

Fitness Star is a rule based expert system built on JESS that is designed to take certain health related parameters from the user and infer certain vital parameters and makes recommendations to the user of the system. The system advices user his workout and food intake plans based on the person's age, gender, weight, and height and resting heart rate. The system is also capable of diagnosing chances of diabetes and other coronary heart conditions based on blood pressure levels, sugar levels and sleep patterns.

#### Features:

- 1. The system can calculate parameters like Body Mass Index and Basal Metabolic Rate of the individual.
- 2. The system calculates Body Fat percentage and Lean Body Mass.
- 3. The system can classify the individual into different obesity fitness levels.
- 4. The system also computes the current daily food intake in calories and recommends to increase or decrease food intake depending on certain vital parameters.
- 5. The system also calculates the number of calories to be burnt per day for individuals who are above the normal weight.
- 6. The system also gives you of a rough estimate of number of days it will take to reach the desired fitness weight, subject to the user pursues the recommendations of the system.
- 7. The system is also capable to identifying stress based on hypertension symptoms and sleep patterns.
- 8. Based on obesity indices and stress levels the system also estimates risk of Coronary heart diseases.
- 9. The system also infers chances of diabetes.
- 10. In case the system finds users with extreme symptoms, the system recommends the user to consult a physician immediately.
- 11. The system exits by displaying a different motivational quote in the end.

# Rules and descriptions:

#	Rule Name	Description	References and
1	calculateBMI	Calculates the BMI based	detailed explanations BMI Wiki
_	our	on height and weight	DIVIT VVIKI
2	checkBMIRating	Assigns a BMI rating.	BMI Wiki
3	assignFatLvlDesc	Assigns an additional	BMI Wiki
		obesity comment.	
4	calcBodyFatPercent	Calculates the percentage	Calculating Fat
		of fat in the body	<u>Percentage</u>
5	calcLeanBodyMass	Calculate the lean body	<u>Lean Body Mass Wiki</u>
		mass of the individual	
6	weightLossGainPlan	Plan and estimate	Using the BMI formula
		parameters like weight to	with the desired BMI.
		be gained or lost.	
7	calculateBMR	Calculates Basal	BMR formula
		Metabolic Rate of an	
		individual.	
8	calcCalIntakeToMaint	Calculate the current	<u>Harris Benedict</u>
		daily average calorie	<u>Equation</u>
		intake (Daily calorie	
		needs to maintain current	
		weight) based on Harris	
	and a Call Table by Down Ober	Benedict Formula	
9	calcCalIntakeBurnQty	Calculate the number of	U.S. Army fitness
		calories to intake so as to	manuals recommends
		lose or gain weight.	reducing intake by 15%
10	createWorkoutPlan	Calculates the total	to assist weight loss.
10	CleateworkoutFlan	number of calories to be	Assumption the person will run at 6
		burnt to reach target	mph(moderate) or 10
		weight and calculate the	mph(intense) and will
		number of days to reach	only run up to a
		out the desired weight,	maximum of 90
		subject to the user	minutes
		pursuing the plan	advised by American
		vigorously.	College of Sports
		,	Medicine
11	calcCalIntakeGainQty	Calculate the new intake	U.S. Army fitness
		calories for under-weight	manuals recommends
		individuals to gain the	increasing the intake by
		desired weight.	500 calories daily to
			assist weight gain.
12	computeTrgtDaysWgtGain	Calculates the number of	
		days to achieve the	
		desired weight for	
		underweight individuals.	

12	assessDiabeticCondition	Dia ana a a a a a a dha a dha a a a a	Consentational
13	assessDiabeticCondition	Diagnoses whether an	Sugar levels and
		individual has any	diabetes
	checkBldPressLvls	symptoms of diabetes.	
14	CheckBldFressLvIs	Checks and assigns	<u>Understanding blood</u>
		suitable hypertension	pressure readings
	1 10: - 1	levels to the individual.	
15	checkStrsLvls	Checks if the individual is	
		stressed based on hyper	
		tension levels and sleep	
	1 1 0 0 0 1	patterns.	
16	chckCHDRisk	Check for any coronary	
		heart disease symptoms	
		based on stress levels and	
		hyper tension	
17	adjustForObesity	Take into account obesity	
		and magnify it if	
		necessary.	
18	calcWaterQty	Calculate the ideal	Calculating amount of
		amount of water intake	water intake every day
		depending on personal	
		parameters for a healthy	
		life style.	
19	findTargetHR	Compute the target heart	<u>Calculating your target</u>
		rate to be achieved	<u>heart rate</u>
		during the workout	
		sessions based on	
		person's age, average	
		heart rate and workout	
		intensity.	
20	chckConsultingReq	In case of extreme	Emergency scenarios
		scenarios recommend for	
		immediate physician	
		consulting.	
The	below rules reside in the <b>output.clp</b> file and	are only used to direct outpu	t to terminal has no
proc	essing significance.		
21	dispOutputWelcome	Displays a welcome banner and personal	
		information.	
22	disImmediateDiagnosisDiab	Displays any diagnosis related to diabetes.	
23	disImmediateDiagnosisHeartCondn	Displays any diagnosis related to coronary heart	
		conditions.	
24	fitnessOutput	Displays all outputs related to fitness parameters,	
L		recommendations for calorie intake etc.	
25	workoutOutput	Displays all workout parameters in case the person	
		is over-weight.	

## Usage Manual:

#### Instructions:

Copy the files health.clp, input.clp and output.clp to the BIN folder under the JESS directory.

Open JESS and execute the below commands:

```
(batch health.clp)
(batch input.clp)
(batch output.clp)
```

In case the grader wants to change inputs and test new inputs please only make changes to the input.clp file. Also, make sure there is only one active case in each run i.e. the input.clp file will have only one assert and one run statement at a time.

### Sample runs:

#### Run #1

#### The contents of input.clp file:

(assert (person(personName Ironman) (sex M) (activityType light) (sugarLevel
100.6) (workoutType moderate) (age 24) (height 2.0) (weight
95) (bloodPressureSystolic 118) (bloodPressureDiastolic 81) (currentSleepHours
2) (avgHR 125)))

(run)

```
D:\less7ip2\bino:jess.bat

Jess, the Rule Engine for the Java Platform
Copyright (C) 2008 Sandia Corporation
Jess Version 7.1p2 11/5/2008

Jess> (batch health.clp)

TRUE
Jess> (batch input.clp)

Jess> (batch input.clp)

Jess> (batch output.clp)

Jess> (batch output.clp)

Jess> (batch output.clp)

Jess> (batch output.clp)

Mello there. Ironman!

Below is the information you entered.
Your gender is You regular life style is a light lifestyle

Based on your hight and weight your Body Mass Index is 23.75

Based on your hight and weight your Body Mass Index is 23.75

Based on your MIY you have been categorized as Normal Range
Your body fat percentage is 17.820000000000004

Based on your sugar level I diagnosed You have High chances of diabetes
I recommend you to consult your physician in order for proper medication if required.

Based on your blood pressure levels, sleep patterns; obesity levels and stress levels I have diagnosed your coronary heart condition as less
I recommend you to consult your physician in order for proper medication if required.

I think you are already fit!

Hope to see you again.
Thank you for using the Fitness Star!
Keep calm and Keep working out!

O:\Jess71p2\bino

Jess 11 | August 11 | August 12 | August 12 |

Jess 2000 | August 12 | August 12 |

Jess 2000 | August 12 |

Jess 2000 | August 12 |

Jess 21 | August 12 |

Jess 22 | August 12 |

Jess 21 | August 12 |

Jess 22 | August 12 |

Jess 22 | August 12 |

Jess 22 | August 12 |

Jess 23 | August 12 |

Jess 24 | August 12 |

Jess 25 | August
```

#### Run #2

#### The contents of input.clp file:

```
(assert (person(personName Catwoman) (sex F) (activityType
sedentary) (sugarLevel 100.6) (workoutType moderate) (age 24) (height
1.6256) (weight 40) (bloodPressureSystolic 125) (bloodPressureDiastolic
90) (currentSleepHours 5) (avgHR 125)))
(run)
```

```
Divisor/Ip2(binn-yess.bat

Jess, the Rule Engine for the Java Platform
Copyright (C) 2008 Sandia Corporation
Jess Version 7.1p2 11/5/2008

Jess> (batch health.clp)

THE
Jess> (batch input.clp)

Jess> (batch input.clp)

Jess> (batch output.clp)

Jess (batch outp
```

#### **Run #3**

#### The contents of input.clp file:

(assert (person(personName Batman) (sex M) (activityType intense) (sugarLevel 120.0) (workoutType moderate) (age 24) (height 1.6256) (weight 90) (bloodPressureSystolic 140) (bloodPressureDiastolic 101) (currentSleepHours 9) (avgHR 100))) (run)

```
Divises/Tap2\binojess.bat

Jess, the Rule Engine for the Java Platform
Copyright (C) 2008 Sandia Corporation
Jess Version 7.1p2 11/5/2008

Hello there. Batmani

Below is the information you entered.

Your good and the set of the
```

#### Test cases:

#### Please use only one test case at a time in the input.clp file

```
(assert (person(personName Robin) (sex M) (activityType sedentary) (sugarLevel
65.6) (workoutType moderate) (age 24) (height 1.7272) (weight
76) (bloodPressureSystolic 125) (bloodPressureDiastolic 85) (currentSleepHours
8) (avgHR 125)))
(assert (person(personName Batman) (sex M) (activityType intense) (sugarLevel
120.0) (workoutType moderate) (age 24) (height 1.6256) (weight
90) (bloodPressureSystolic 140) (bloodPressureDiastolic 101) (currentSleepHours
9) (avgHR 100)))
(assert (person(personName Catwoman) (sex F) (activityType
sedentary) (sugarLevel 100.6) (workoutType moderate) (age 24) (height
1.6256) (weight 40) (bloodPressureSystolic 125) (bloodPressureDiastolic
90) (currentSleepHours 5) (avgHR 125)))
(assert (person(personName Ironman) (sex M) (activityType light) (sugarLevel
100.6) (workoutType moderate) (age 24) (height 2.0) (weight
95) (bloodPressureSystolic 118) (bloodPressureDiastolic 81) (currentSleepHours
2) (avgHR 125)))
```